


Product Name: Recombinant Human SMAD6 Protein
Catalog #: PEH7394



Summary

Name	SMAD6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43541
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RRLWRSRVVP DREEGGSGGG GGGDEDGSLG SRAEPAPRAR EGGGCGRSEV RPVAPRRPRD AVGQRGAQGA GRRRRRAGGPP

image



Background

Alternative Names	Mothers against decapentaplegic homolog 6 (MAD homolog 6) (Mothers against DPP homolog 6) (SMAD family member 6) (SMAD 6) (Smad6) (hSMAD6)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use