


Product Name: Recombinant Human SNAI2 Protein
Catalog #: PEH7396



Summary

Name	SNAI2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43623
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PSDTSSKDHS GSESPISDEE ERLQSKLSDP HAIEAEKFQC NLCNKTYSTF SGLAKHKQLH CDAQSRKSFS CKYCDKEYVS

image



Background

Alternative Names Zinc finger protein SNAI2 (Neural crest transcription factor Slug) (Protein snail homolog 2)

Note

For Research Use Only , Not for Diagnostic Use