


Product Name: Recombinant Human TCF21 Protein
Catalog #: PEH7397



Summary

Name	TCF21_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43680
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FVTSNESTEE SSNCENGSPQ KGRGGLGKRR KAPTKKSPLS GVSQEGKQVQ RNAANARERA RMRVLSKAFS RLKTTLPWVP PDTKLSKLDL LRLASSYIAH

image



Background

Alternative Names Transcription factor 21 (TCF-21) (Capsulin) (Class A basic helix-loop-helix protein 23) (bHLHa23) (Epicardin) (Podocyte-expressed 1) (Pod-1)

Note

For Research Use Only , Not for Diagnostic Use