


Product Name: Recombinant Human AP1G1 Protein
Catalog #: PEH7398



Summary

Name	AP1G1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43747
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IKRRAMELSF ALVNGNNIRG MMKELLYFLD SCEPEFKADC ASGIFLAAEK YAPSKRWHID TIMRVLTTAG SYVRDDAVPN LIQLITNSVE MHAYTVQRLY KAILGDYSQQ PLVQVAAWCI GEYGDLLVSG QCEEEPIQV TEDEVLDILE SVLISNMSTS VTRGYALTAI MKLSTRFTCT VNRIKKVVS I YGSSIDVELQ QRAVEYNALF KKYDHMRSAL LERMPVMEKV TTNGPTEIVQ TNGETEPAPL ETKPPPSGPQ PTSQANDLLD LLGGNDITPV IPTAPTSKPS SAGGELDLL GDINLTGAPA

image



Background

Alternative Names	AP-1 complex subunit gamma-1 (Adaptor protein complex AP-1 subunit gamma-1) (Adaptor-related protein complex 1 subunit gamma-1) (Clathrin
--------------------------	---

Product Name: Recombinant Human AP1G1 Protein
Catalog #: PEH7398



assembly protein complex 1 gamma-1 large chain) (Gamma1-adaptin) (Golgi
adaptor HA1/AP1 adaptin subunit gamma-1)

Note

For Research Use Only , Not for Diagnostic Use