

Product Name: Recombinant Human NRDC Protein
Catalog #: PEH7401



Summary

Name	NRDC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43847
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TQATKYNSEV VDKKIEEFLS SFEEKIENLT EEFNTQVTA LIKLKECEDT HLGEEVDRNW NEVVTQQYLF DRLAHEIEAL KSFSKSDLVN WFKAHRGPGS KMLSVHVVG YGKYEEEDGT PSSSDSNSSC EVMQLTYLPT SPLADCIIP

image

SDS

Background

Alternative Names	Nardilysin (EC 3.4.24.61) (N-arginine dibasic convertase) (NRD convertase) (NRD-C) (Nardilysin convertase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use