


Product Name: Recombinant Human AHCYL1 Protein
Catalog #: PEH7402



Summary

Name	SAHH2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43865
Host	E.coli
Species	Human
Molecular Mass	60 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSMPDAMPLP GVGEELKQAK EIEDAEKYSF MATVTKAPKK QIQFADDMQE FTKFPTKTGR RLSRSISQS STDSYSSAAS YTDSSDDEVS PREKQQTNSK GSSNFCVKNI KQAEFGRREI EIAEQDMSAL ISLRKRAQGE KPLAGAKIVG CTHITAQTAV LIETLCALGA QCRWSACNIY STQNEVAAAL AEAGVAVFAW KGESEDDFWW CIDRCVNMDG WQANMILDDG GDLTHWVYKK YPNVFKKIRG IVEESVTGVH RLYQLSKAGK LCVPMNVND SVTKQKFDNL YCCRESILDG LKRTTDMVFG GKQVVVCGYG EVGKGCCAAL KALGAIVYIT EIDPICALQA CMDGFRVVKL NEVIRQVDVV ITCTGNKNVV TREHLDRMKN SCIVCNMGHS NTEIDVTSR TPETWVERV SQVDHVIWPD GKRVLVLAEG RLLNLSCSTV PTFVLSITAT TQALALIELY NAPEGRYKQD VYLLPKKMDE YVASLHLPSP DAHLTELTDD QAKYLGLNKN GPFKPNYYRY

image



Background

Product Name: Recombinant Human AHCYL1 Protein
Catalog #: PEH7402



Alternative Names

S-adenosylhomocysteine hydrolase-like protein 1 (DC-expressed AHCY-like molecule) (IP(3)Rs binding protein released with IP(3)) (IRBIT) (Putative adenosylhomocysteinase 2) (S-adenosyl-L-homocysteine hydrolase 2) (AdoHcyase 2)

Note

For Research Use Only , Not for Diagnostic Use