


Product Name: Recombinant Human KALRN Protein
Catalog #: PEH7403



Summary

Name	KALRN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60229
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LWYLRLLRL DSGSFRNDGL KASDVLPILK EKVAFVSGGR DKRGGPILTF PARSNHDRIR QEDLRKLVTY LASVPSEDVC KRGFTVIIDM RGSKWDLIKP LLKTLQEAFP AEIHVALIHK PDNFWQKQKT NFGSSKFIFE TSMVSVEGLT KLVDPSQLTE EFDGSLDYNH EEWIELRLSL EEFNSAVHL LSRLEDLQEM LARKEFPVDV EGSRRLLIDEH

image



Background

Alternative Names	Kalirin (EC 2.7.11.1) (Huntingtin-associated protein-interacting protein) (Protein Duo) (Serine/threonine-protein kinase with Dbl- and pleckstrin homology domain)
--------------------------	--

Product Name: Recombinant Human KALRN Protein
Catalog #: PEH7403



Note

For Research Use Only , Not for Diagnostic Use