


Product Name: Recombinant Human ZEB2 Protein
Catalog #: PEH7404



Summary

Name	ZEB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60315
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FQGELLDGAQ DYMSGGLDDMT DSDSCLSRKK IKKTESGMYA CDLCDKTFQK SSSLRHKYE HTGKRPHQCQ ICKKAFKHKH HLIHSRLHS GEKPYQCDKC GKRFSHSGSY SQHMNHRYSY CKREAEEREA AEREAREKGH LEPTLLMNR AYLQSITPQG YSDSEERESM PRDGESEKEH EKEGEDGYGK LGRQDGDDEEF EEEEEESENK SMDTDPETIR DEEETGDHSM DDSSEDGKME TKSDHEEDNM

image



Background

Alternative Names	Zinc finger E-box-binding homeobox 2 (Smad-interacting protein 1) (SMADIP1) (Zinc finger homeobox protein 1b)
--------------------------	---

Note

Product Name: Recombinant Human ZEB2 Protein
Catalog #: PEH7404



For Research Use Only , Not for Diagnostic Use