


Product Name: Recombinant Human SNAP91 Protein
Catalog #: PEH7405



Summary

Name	AP180_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60641
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PEKLLKSMP I LQGQIDALLE FDVHPNELTN GVINA AFMLL FKDLIKLFAC YNDGVINLLE KFFEMKKGQC KDALEIYKRF LTRMTRVSEF LKVAEQVGID KGDIPDLTQA PSSLMETLEQ HLNTLEGKKP GNNEGSGAPS

image



Background

Alternative Names	Clathrin coat assembly protein AP180 (91 kDa synaptosomal-associated protein) (Clathrin coat-associated protein AP180) (Phosphoprotein F1-20)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use