


Product Name: Recombinant Human ATP6V1G1 Protein
Catalog #: PEH7407




Summary

Name	VATG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75348
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MASQSQGIQQ LLQAEKRAAE KVSEARKRKN RRLKQAKEEA QAEIEQYRLQ REKEFKAKEA AALGSRGSCS TEVEKETQEK MTILQTYFRQ NRDEVLDNLL AFVCDIRPEI HENYRING

image



SDS



WB

Background

Alternative Names	V-type proton ATPase subunit G 1 (V-ATPase subunit G 1) (V-ATPase 13 kDa subunit 1) (Vacuolar proton pump subunit G 1) (Vacuolar proton pump subunit M16)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ATP6V1G1 Protein
Catalog #: PEH7407

