

Product Name: Recombinant Human ULK1 Protein
Catalog #: PEH7411



Summary

Name	ULK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75385
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GSSSSSSSTS HLASPPSLGE MQQLQKTLAS PADTAGFLHS SRDSGGGSKDS SCDTDDFVMV PAQFPGDLVA EAPSAKPPPD SLMCSGSSLV ASAGLESHGR TPSPSPPCSS SPSPSGRAGP FSSSRCGASV

image

SDS

Background

Alternative Names	Serine/threonine-protein kinase ULK1 (EC 2.7.11.1) (Autophagy-related protein 1 homolog) (ATG1) (hATG1) (Unc-51-like kinase 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use