


Product Name: Recombinant Human RPS6KA4 Protein
Catalog #: PEH7414



Summary

Name	KS6A4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75676
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AVELRITEAN LTGHEEKVSV ENFELLKVLG TGAYGKVFLV RKAGGHDAGK LYAMKVLKKA ALVQRAKTQE HTRTERSVLE LVRQAPFLVT LHYAFQTDAG LHLILDYVSG GEMFTHLYQR QYFKEAEVRV YGGEIVLAL HLHKLGIYR DLKLENVLLD SEGHIULTDF GLSKEFLTEE KERTFSFCGT IEYMAPEIIR SKTGHGKA VD WWSLGILLFE LLTGASPFTL EGERNTQAEV SRRILKCSPP FPPRIGPVAQ DLLQRLCKD PKKRLGAGPQ GAQEVNRNHPF FQGLDWVALA ARKIPAPFRP QIRSELDVGN FAEFTRLEP VYSPPGSPPP GDPRIFQGY FVAPSILFDH NNAVMTDGL APGAGDRPGR AAVARSAMMQ

image



Background

Alternative Names	Ribosomal protein S6 kinase alpha-4 (S6K-alpha-4) (EC 2.7.11.1) (90 kDa)
--------------------------	--

Product Name: Recombinant Human RPS6KA4 Protein
Catalog #: PEH7414



ribosomal protein S6 kinase 4) (Nuclear mitogen- and stress-activated
protein kinase 2) (Ribosomal protein kinase B) (RSKB)

Note

For Research Use Only , Not for Diagnostic Use