


Summary

Name	STK10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94804
Host	E.coli
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEQMHRFEQ EINAKKKFFD TELENLERQQ KQQVEKMEQD HAVRRREEAR RIRLEQDRDY TRFQEQLKLM KKEVKNEVEK LPRQQRKESM KQKMEEHTQK KQLLDRDFVA KQKEDLELAM KRLTTDNRRE ICDKERECLM KKQELLRDRE AALWEMEEHQ LQERHQLVKQ QLKDQYFLQR HELLRKHEKE REQMQRYNQR MIEQLKVRQQ QEKARLPKIQ RSEGKTRMAM YKKS LHINGG GSAAEQREKI KQFSQQEEKR QKSERLQQQQ KHENQMRDML AQCESNMSEL QQLQNEKCHL LVEHETQKLIK ALDESHNQNL KEWRDKLRPR KKALEEDLNQ KKREQEMFFK LSEEAECNP STPSKAAKFF

image



Background

Alternative Names Serine/threonine-protein kinase 10 (EC 2.7.11.1) (Lymphocyte-oriented

Product Name: Recombinant Human STK10 Protein
Catalog #: PEH7419



kinase)

Note

For Research Use Only , Not for Diagnostic Use