

Summary

Name	SC31A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94979
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RMPSHQGAEQ QQQQHVFIS QVVTEKEFLS RSDQLQQAVQ SQGFINYCQK KIDASQTEFE KNVWSFLKVN FEDDSRGKYL ELLGYRKEDL GKKIALALNK VDGANVALKD SDQVAQSDGE ESPAAEEQLL GEHIKEEKEE SEFLPSSGGT FNISVSGDID GLITQALLTG NFESAVDLCL HDNRMADAI

image

SDS

Background

Alternative Names	Protein transport protein Sec31A (ABP125) (ABP130) (SEC31-like protein 1) (SEC31-related protein A) (Web1-like protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use