


**Product Name: Recombinant Human ATE1 Protein**  
**Catalog #: PEH7422**



## Summary

<b>Name</b>	ATE1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O95260
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	60 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MAFWAGGSPS VVDYFPSEDF YRCGYCKNES GSRNMGMAH SMTVQDYQDL IDRGWRRSGK YVYKPVMNQT CCPQYTIRCR PLQFQPSKSH KKVLLKMLKF LAKGEVPGKS CEDEPMDSTM DDAVAGDFAL INKLDIQCDL KTLSDDIKES LESEGKNSKK EEPQELLQSQ DFVGEKLGSG EPSHSVHVHT VPKPGKGADL SKPPCRKAKE IRKERKRLKL MQQNPAGELE GFQAQGHPPS LFPPKAKSNQ PKSLEDLIFE SLPENASHKL EVRVVRSPP SSQFKATLLE SYQVYKRYQM VIHKNPPDTP TESQFTRFLC SSPLEAETPP NGPDCGYGSF HQQYWLDGKI Iavgvidilp NCVSSVYLY DPDYSFSLG VYSALREIAF TRQLHEKTSQ LSYYYMGFYI HSCPMMKYKG QYRPSDLLCP ETYVWVPIEQ CLPSLENSKY CRFNQDPEAV DEDRSTEPDR LQVFHKRAIM PYGVYKKQQK DPSEEAVALQ YASLVGQKCS ERMILLFRN

## image



## Background

**Product Name: Recombinant Human ATE1 Protein**  
**Catalog #: PEH7422**



---

**Alternative Names**

Arginyl-tRNA--protein transferase 1 (Arginyltransferase 1) (R-transferase 1)  
(EC 2.3.2.8) (Arginine-tRNA--protein transferase 1)

**Note**

For Research Use Only , Not for Diagnostic Use