


Product Name: Recombinant Human TACC2 Protein
Catalog #: PEH7423



Summary

Name	TACC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95359
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQIARAEIIT KEREVSEWKD KYEESRREVM EMRKIVAEYE KTIAQMIEDE QREKSVSHQT VQQLVLEKEQ ALADLNSVEK SLADLFRRYE KMKEVLEGFR KNEEVLRCA QEYLSRVKKE EQRYQALKVH AEEKLDRANA EIAQVRGKAQ QEQAHHQASL RKEQLRVDAL ERTLEQKNKE IEELTKICDE

image



Background

Alternative Names Transforming acidic coiled-coil-containing protein 2 (Anti-Zuai-1) (AZU-1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TACC2 Protein
Catalog #: PEH7423

