


Product Name: Recombinant Human ADCY5 Protein
Catalog #: PEH7426



Summary

Name	ADCY5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95622
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GYAAQKTAAP APRGGPEHRS AWGEADSRAN GYPHAPGGSA RGSTKKPGGA VTPQQQQLRL SRWRSDDDDD PPLSGDDPLA GGFGFSFRSK SAWQERGGDD CGRGSRQRR GAASGGSTRA PPAGGGGGSA AAAASAGGTE VRPRSVEVGL EERRGKGRAA DELEAGAVEG GEGSGDGGSS

image



Background

Alternative Names	Adenylate cyclase type 5 (EC 4.6.1.1) (ATP pyrophosphate-lyase 5) (Adenylate cyclase type V) (Adenylyl cyclase 5) (AC5)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ADCY5 Protein
Catalog #: PEH7426

