


Product Name: Recombinant Human EML2 Protein
Catalog #: PEH7432



Summary

Name	EMAL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	O95834
Host	E.coli
Species	Human
Molecular Mass	68 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YLGHNDDIKC LAIHPDMVTI ATGQVAGTTK EGKPLPPHVR IWDSVSLSTL HVLGLGVFDR AVCCVGFSKS NGGNLLCAVD ESNDHMLSVW DWAKETKVVD VKCSNEAVLV ATFHPTDPTV LITCGKSHIY FWTLEGGSLK KRQGLFEKHE KPKYVLCVTF LEGGDVVTGD SGGNLYVWGK GGNRITQAVL GAHDGGVFGL CALRDGTLVS GGGRDRRVVL WGS DYSLQVE VEPEDFGPV RTVAEGHGDT LYVGTTNSI LQGSVHTGFS LLVQGHVEEL WGLATHPSRA QFVTCGQDKL VHLWSSDSHQ PLWSRIEDP ARSAGFHPSG SVLAVGTVTG RWLLLDTETH DLVAIHTDGN EQISVVSFSP DGAYLAVGSH DNLVYVYTVD QGGRKVSRLG KCSGHSSFIT HLDWAQDSSC FVTNSGDYEI LYWDPATCKQ ITSADAVRNM EWATATCVLG FGVFGIWSEG ADGTDINAVA RSHDGKLLAS

image



Background

Product Name: Recombinant Human EML2 Protein
Catalog #: PEH7432



Alternative Names

Echinoderm microtubule-associated protein-like 2 (EMAP-2) (HuEMAP-2)

Note

For Research Use Only , Not for Diagnostic Use