


Summary

Name	BPNT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95861
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MASSNTVLMR LVASAYSIAQ KAGMIVRRVI AEGDLGIVEK TCATDLQTKA DRLAQMSICS SLARKFPKLT IIGEEDLPSE EVDQELIEDS QWEEILKQPC PSQYSAIKEE DLVWVVDPLD GTKEYTEGLL DNVTVLIGIA YEGKAIAGVI NQPYYNYEAG PDAVLGRTIW GVLGLGAFGF QLKEVPAGKH IITTRSHSN KLVTDCAAM NPDAVLRVGG AGNKIIQLIE GKASAYVFAS PGCKKWDTCA PEVILHAVGG KLTDIHGNVL QYHKDVKHMN SAGVLATLRN YDYYASRVPE SIKNALVP

image



Background

Alternative Names	3'(2'),5'-bisphosphate nucleotidase 1 (EC 3.1.3.7) (Bisphosphate 3'-nucleotidase 1) (PAP-inositol 1,4-phosphatase) (PIP)
--------------------------	--

Product Name: Recombinant Human BPNT1 Protein
Catalog #: PEH7433



Note

For Research Use Only , Not for Diagnostic Use