


Product Name: Recombinant Human KRT1 Protein
Catalog #: PEH7440



Summary

Name	K2C1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P04264
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	REQIKSLNNQ FASFDKVRF LEQQNQLVLT KWELLQQVDT STRTHNLEPY FESFINNLR RVDQLKSDQS RLDSELKNMQ DMVEDYR NKY EDEINKRTNA ENEFVTIKKD VDGAYMTKVD LQAKLDNLQQ EIDFLTALYQ AELSQQMTQI SETNVILSMD NNRSLLDSI IAEVKAQYED IAQKSKAEAE SLYQSKYEEL QITAGRHGDS VRNSKIEISE LNRVIQLRS EIDNVKKQIS NLQQSISDAE QRGENALKDA KNKLNLEDA LQQAKEDLAR LLRDYQELMN TKLALDLEIA TYRTLLEGEE SRMSGECAPN VSVSVSTSH

image



Background

Alternative Names Keratin, type II cytoskeletal 1 (67 kDa cytokeratin) (Cytokeratin-1) (CK-1) (Hair alpha protein) (Keratin-1) (K1) (Type-II keratin Kb1)

Product Name: Recombinant Human KRT1 Protein
Catalog #: PEH7440



Note

For Research Use Only , Not for Diagnostic Use