


Product Name: Recombinant Human HMGN2 Protein
Catalog #: PEH7446



Summary

Name	HMGN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P05204
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MPKRKAEGDA KGDKAKVKDE PQRRLSA KPAPPKPEPK PKKAPAKKGE KVPKGKKGKA DAGKEGNNPA ENGDAKTDQA QKAEGAGDAK

image



Background

Alternative Names Non-histone chromosomal protein HMG-17 (High mobility group nucleosome-binding domain-containing protein 2)

Note

For Research Use Only , Not for Diagnostic Use