


Product Name: Recombinant Human TFAP2A Protein
Catalog #: PEH7447



Summary

Name	AP2A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P05549
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDLLHGPHAL SSGLGDSIH SLPHAIEEVP HVEDPGINIP DQTVIKKGPV SLSKSNSNAV SAIPINKDNL FGGVVPNEV FCSVPGRSL LSSTSKYKVT VAEVQRRLLSP PECLNASLLG GVLRRASKN GGRSLREKLD KIGLNLPAQR RKAANVTLLT SLVEGEAVHL ARDFGYVCET EFPKAVAEF LNRQHSDPNE QVTRKNMLLA TKQICKEFTD LLAQDRSPLG NSRPNPILEP GIQSCLTHTN LISHGFGSPA VCAAVTALQN YLTEALKAMD KMYLSNNPNS HTDNNAKSSD

image



Background

Alternative Names	Transcription factor AP-2-alpha (AP2-alpha) (AP-2 transcription factor) (Activating enhancer-binding protein 2-alpha) (Activator protein 2) (AP-2)
--------------------------	--

Product Name: Recombinant Human TFAP2A Protein
Catalog #: PEH7447



Note

For Research Use Only , Not for Diagnostic Use