


Product Name: Recombinant Human FH Protein
Catalog #: PEH7451



Summary

Name	FUMH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P07954
Host	E.coli
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASQNSFRIEY DTFGELKVPN DKYYGAQTVR STMNFKIGGV TERMPTPIK AFGILKRAAA EVNQDYGLDP KIANAIMKAA DEVAEGKLND HFPLVVWQTG SGTQTNMNVN EVISNRAIEM LGGELGSKIP VHPNDHVNKS QSSNDFPTA MHIAAAIEVH EVLLPGLQKL HDALDAKSKE FAQIIKIGRT HTQDAVPLTL GQEFSGYVQQ VKYAMTRIKA AMPRIYELAA GGTAVGTGLN TRIGFAEKVA AKVAALTGLP FVTAPNKFEA LAAHDALVEL SGAMNTTACS LMKIANDIRF LGSGPRSLG ELILPENEPG SSIMPGKVNP TQCEAMTMVA AQVMGNHVAV TVGGSNGHFE LNVFKPMMIK NVLHSARLLG DASVSFTENC VVGIQANTER INKLMNESLM LVTALNPHIG YDKAAKIAKT AHKNGSTLKE TAIELGYLTA EQFDEWVKPK DMLGPK

image



Background

Product Name: Recombinant Human FH Protein
Catalog #: PEH7451



Alternative Names Fumarate hydratase, mitochondrial (Fumarase) (HsFH) (EC 4.2.1.2)

Note

For Research Use Only , Not for Diagnostic Use