


Product Name: Recombinant Human VIL1 Protein
Catalog #: PEH7454



Summary

Name	VILI_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P09327
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	THTVGSVAKV EQVKFDATSM HVKPQVAAQQ K MVDDGSGEV QVWRIENLEL VPVDSKWLGH FYGGDCYLLL YTYLIGEKQH YLLYVWQGSQ ASQDEITASA YQAVILDQKY NGEPVQIRVP MGKEPPHLS IFKGRMVVYQ GGTSRTNNLE TGPSTRLFQV QGTGANNTKA FEVPARANFL NSNDVFVLKT QSCCYLWCGK GCSGDEREMA K MVADTISR T EKQVVVEGQE PANFWMALGG KAPYANTKRL QEENLVITPR LFEC SNKTGR FLATEIPDFN QDDLEEDDFV LLDVWDQVFF WIGKHANEEE K KAAATTAQE YLKTHPSGRD PETPIIVVKQ GHEPPTFTGW FLAWDPFKWS

image



Background

Alternative Names Villin-1

Product Name: Recombinant Human VIL1 Protein
Catalog #: PEH7454



Note

For Research Use Only , Not for Diagnostic Use