

**Product Name: Recombinant Human IFIT1 Protein**  
**Catalog #: PEH7456**



## Summary

<b>Name</b>	IFIT1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P09914
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	56 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MSTNGDDHQV KDSLEQLRCH FTWELSIDDD EMPDLENRVL DQIEFLDTKY SVGIHNLLAY VKHLKGQNEE ALKSLKEAEN LMQEEHDNQA NVRSVLTWGN FAWMYYHMGR LAEAQTYLDK VENICKKLSN PFRYRMECPE IDCEEGWALL KCGGKNYERA KACFEKVLEV DPENPESSAG YAISAYRLDG FKLATKNHKP FLLPLRQAV RLNPDNGYIK VLLALKLQDE GQEAEGEKYI EEALANMSSQ TYVFRYAAKF YRRKGSVDKA LELLKALQE TPTSLLHHQ IGLCYKAQMI QIKEATKGQP RGQNREKLDK MIRSAIFHFE SAVEKKPTFE VAHLDLARMY IEAGNHRKAE ENFQKLLCMK PVVEETMQDI HFHYGRFQEF QKKSDVNAII HYLKAIKIEQ ASLTRDKSIN SLKKLVLRLK RRKALDLESL SLLGFVYKLE GNMNEALEYY ERALRLAADF ENSVRQGP

## image



## Background

**Product Name: Recombinant Human IFIT1 Protein**  
**Catalog #: PEH7456**



---

**Alternative Names**

Interferon-induced protein with tetratricopeptide repeats 1 (IFIT-1)  
(Interferon-induced 56 kDa protein) (IFI-56K) (P56)

**Note**

For Research Use Only , Not for Diagnostic Use