


Product Name: Recombinant Human H1-4 Protein
Catalog #: PEH7459



Summary

Name	H14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P10412
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APAPAEKTPV KKKARKSAGA AKRKASGPPV SELITKAVAA SKERSGVSLA ALKKALAAAG YDVEKNSRI KLGLKSLVSK GTLVQTKGTG ASGSFKLNKK AASGEAKPKA KKAGAAKAKK PAGAACKPKK ATGAATPKKS AKKTPKKAKK PAAAAGAKKA

image



Background

Alternative Names Histone H1.4 (Histone H1b) (Histone H1s-4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human H1-4 Protein
Catalog #: PEH7459

