



Product Name: Recombinant Human EPB41 Protein
Catalog #: PEH7462



Summary

Name	EPB41_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P11171
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HCKVSLDDT VYECVVEKHA KGQDLLKRV EHLNLEEDY FGLAIWDNAT SKTWLDSAKE IKKQVRGVPW NFTFNVKFYP PDPAQLTEDI TRYLLCLQLR QDIVAGRLPC SFATLALLGS YTIQSELGDY DPELHGVDYV SDFKLAPNQT KELEEKVMEL HKSYSRSMTPA QADLEFLNA KKLSMYGVDL HKAKDLEGVD IILGVCSSGL LVIKDKLRIN RFPWPKVLKI SYKRSSFFIK IRPGEQEQYE STIGFKLPSY RAAKKLWKVC VEHTFFRLT

image



Background

Alternative Names Protein 4.1 (P4.1) (4.1R) (Band 4.1) (EPB4.1) (Erythrocyte membrane protein band 4.1)

Product Name: Recombinant Human EPB41 Protein
Catalog #: PEH7462



Note

For Research Use Only , Not for Diagnostic Use