


Product Name: Recombinant Human PABPC1 Protein
Catalog #: PEH7468



Summary

Name	PABP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P11940
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASLYVGDLHP DVTEAMLYEK FSPAGPILSI RVC RDMITRR SLGYAYVNFQ QPADAERALD TMNFDVIK GK PVRIMWSQRD PSLRKSGVGN IFIKNLDKSI DNKALYDTFS AFGNILSCKV VCDENGSKGY GFVHFETQEA AERAIEKMNG MLLNDRKV FV GRFKSRKERE AELGARAKEF

image



Background

Alternative Names Polyadenylate-binding protein 1 (PABP-1) (Poly(A)-binding protein 1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PABPC1 Protein
Catalog #: PEH7468

