

Product Name: Recombinant Human PDIA4 Protein
Catalog #: PEH7471



Summary

Name	PDIA4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P13667
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSNRENAIED EEEEEEDDD EEEDDLEVKE ENGLVVLNDA NFDNFVADKD TVLLEFYAPW CGHCKQFAPE YEKIANILKD KDPPIPVAKI DATSASVLAS RFDVSGYPTI KILKKGQAVD YEGSRTQEEI VAKVREVSQP DWTPPPEVTL VLTKENFDEV VNDADIILVE FYAPWCGHCK KLAPEYEKAA KELSKRSPPI PLAKVDATEAE TDLAKRFDVS GYPTLKIFRK GRPYDYNGPR EKYGIVDYMI EQSGPPSKEI LTLKQVQEFL KDGDVVIIG VFKGESDPAY QQYQDAANNL REDYKFHHTF STEIAKFLKV SQGQLVVMQP

image

SDS

Background

Alternative Names	Protein disulfide-isomerase A4 (EC 5.3.4.1) (Endoplasmic reticulum resident protein 70) (ER protein 70) (ERp70) (Endoplasmic reticulum resident protein 72) (ER protein 72) (ERp-72) (ERp72)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use