


Product Name: Recombinant Human P4HA1 Protein
Catalog #: PEH7472



Summary

| | |
|------------------------------|--|
| Name | P4HA1_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P13674 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 32 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | LIHTEKDLVT SLKDYIKAE DKLEQIKKWA EKLDRLTSTA TKDPEGFVGH PVNAFKLMKR LNTEWSELEN LVLKDMSDGF ISNLTIQRQY FPNDEDQVGA AKALLRLQDT YNLDTDTISK GNLPGVKHK FLTAEDCFEL GKVAYTEADY YHTELWMEQA LRQLDEGEIS TIDKVSVDY LSYAVYQQGD LDKALLTKK LLELDPEHQR ANGNLYFEY IMAKEKDVNK SASDDQSDQK TTPKKKGAV DYLPERQKYE |

image



Background

Alternative Names Prolyl 4-hydroxylase subunit alpha-1 (4-PH alpha-1) (EC 1.14.11.2)
(Procollagen-proline,2-oxoglutarate-4-dioxygenase subunit alpha-1)

Product Name: Recombinant Human P4HA1 Protein
Catalog #: PEH7472



Note

For Research Use Only , Not for Diagnostic Use