


Product Name: Recombinant Human APEH Protein
Catalog #: PEH7474



Summary

Name	ACPH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P13798
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VDTQVGTVTS LTAGGSGGSW KLLTIDQDLM VAQFSTPSLP PTLKVGFLPS AGKEQSVLWV SLEEAETIPD IHWGIRVLQP PPEQENVQYA GLDFEAILLQ PGSPDKTQV PMVVMPPHGGP HSSFVTAWML FPAMLCKMGF AVLLVNYRGS TGFGQDSILS LPGNVGHQDV KDVQFAVEQV LQEEHFDASH VALMGGSHGG FISCHLIGQY PETYRACVAR NPVINIASML GSTDIPDWCV VEAGFPFSSD CLPDLVWAE MLDKSPIRYI PQVKTPLLLM LGQEDRRVPF KQGMYYRAL KTRNVPVRL LYPKSTHALS EVEVESDSFM

image



Background

Alternative Names	Acylamino-acid-releasing enzyme (AARE) (EC 3.4.19.1) (Acyl-peptide hydrolase) (APH) (Acylaminoacyl-peptidase) (Oxidized protein hydrolase)
--------------------------	--

Product Name: Recombinant Human APEH Protein
Catalog #: PEH7474



(OPH)

Note

For Research Use Only , Not for Diagnostic Use