


Product Name: Recombinant Human JUP Protein
Catalog #: PEH7476



Summary

Name	PLAK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P14923
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HLTSRHPEAE MAQNSVRLNY GIPAIVKLLN QPNQWPLVKA TIGLIRNLAL CPANHAPLQE AAVIPRLVQL LVKAHQDAQR HVAAGTQQPY TDGVRMEEIV EGCTGALHIL ARDPMNRMEI FRLNTIPLFV QLLYSSVENI QRVAAGVLCE LAQDKEAADA IDAEGASAPL MELLHSRNEG TATYAAAVLF RISEDKNPDY RKRVSVELTN SLFKHDPAAW EAAQSMIPIN EPYGDDMDAT YRPMYSSDVP LDPLEMHMDM DGDYPIDTYS DGLRPPYPTA

image



Background

Alternative Names	Junction plakoglobin (Catenin gamma) (Desmoplakin III) (Desmoplakin-3)
--------------------------	--

Product Name: Recombinant Human JUP Protein
Catalog #: PEH7476



Note

For Research Use Only , Not for Diagnostic Use