


Product Name: Recombinant Human TCF3 Protein
Catalog #: PEH7479



Summary

Name	TFE2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P15923
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TAAASEIKRE EKEDEENTSA ADHSEEEKKE LKAPRARTSP DEDEDDLPP EQKAEREKER RVANNARERL RVRDINEAFK ELGRMCQLHL NSEKPQTKLL ILHQAVSVIL NLEQQVRERN LNPKAACKLR REEEKVSGVV

image



Background

Alternative Names	Transcription factor E2-alpha (Class B basic helix-loop-helix protein 21) (bHLHb21) (Immunoglobulin enhancer-binding factor E12/E47) (Immunoglobulin transcription factor 1) (Kappa-E2-binding factor) (Transcription factor 3) (TCF-3) (Transcription factor ITF-1)
--------------------------	--

Note

Product Name: Recombinant Human TCF3 Protein
Catalog #: PEH7479



For Research Use Only , Not for Diagnostic Use