


Product Name: Recombinant Human DSP Protein
Catalog #: PEH7480



Summary

Name	DESP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P15924
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QGIRQPVTVT ELVDSGILRP STVNELESGQ ISYDEVGERI KDFLQGSSCI AGIYNETTKQ KLGIEYAMKI GLVRPGTALE LLEAQAATGF IVDPVSNLRL PVEEAYKRGL VGIEFKEKLL SAERAVTGYN DPETGNIISL FQAMNKELIE KGHGIRLLEA QIATGGIIDP KESHRLPVDI AYKRGYFNEE LSEILSDPSD DTKGFFDPNT EENLTYLQLK ERCIKDEETG LCLLPLKEKK

image



Background

Alternative Names Desmoplakin (DP) (250/210 kDa paraneoplastic pemphigus antigen)

Note

Product Name: Recombinant Human DSP Protein
Catalog #: PEH7480



For Research Use Only , Not for Diagnostic Use