


Product Name: Recombinant Human ALOX15 Protein
Catalog #: PEH7481



Summary

Name	LOX15_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P16050
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GASLYAGSNN QVQLWLVGQH GEALGKRLW PARGKETLK VEVPEYLGPL LFVKLRKRHL LKDDAWFCNW ISVQGGPGAGD EVRFPCYRWV EGNGVLSLPE GTGRTVGEDP QGLFQKHREE ELEERRKLYR WGNWKDGLIL NMAGAKLYDL PVDERFLEDK RVD FEVSLAK GLADLAIKDS LNVLTWCWDL DDFNRIFWCG QSKLAERVRD SWKEDALFGY

image



SDS



WB

Background

Alternative Names	Polyunsaturated fatty acid lipoyxygenase ALOX15 (12/15-lipoyxygenase) (Arachidonate 12-lipoyxygenase, leukocyte-type) (12-LOX) (EC 1.13.11.31) (Arachidonate 15-lipoyxygenase) (15-LOX) (15-LOX-1) (EC 1.13.11.33) (Arachidonate omega-6 lipoyxygenase) (Hepoxilin A3 synthase Alox15) (EC
--------------------------	--

Product Name: Recombinant Human ALOX15 Protein
Catalog #: PEH7481



1.13.11.-) (Linoleate 13S-lipoxygenase) (EC 1.13.11.12)

Note

For Research Use Only , Not for Diagnostic Use