


Product Name: Recombinant Human FER Protein
Catalog #: PEH7483



Summary

Name	FER_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P16591
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HEAVLKLQDW ELRLLETVKK FMALRIKSDK EYASTLQNLC NQVDKESTVQ MNYVSNVSKS WLLMIQQTEQ LSRIMKTHAE DLNSGPLHRL TMMIKDKQQV KKS YIGVHQQ IEAEMIKVTK TELEKLKCSY RQLIKEMNSA KEKYKEALAK GKETEKAKER YDKATMKLHM LHNQYVLALK GAQLHQNQYY DITLPLLLDS LQKMQEEMIK ALKGIFDEYS QITSLVTEEI VNVHKEIQMS

image



Background

Alternative Names	Tyrosine-protein kinase Fer (EC 2.7.10.2) (Feline encephalitis virus-related kinase FER) (Fujinami poultry sarcoma/Feline sarcoma-related protein Fer) (Proto-oncogene c-Fer) (Tyrosine kinase 3) (p94-Fer)
--------------------------	---

Product Name: Recombinant Human FER Protein
Catalog #: PEH7483



Note

For Research Use Only , Not for Diagnostic Use