


Product Name: Recombinant Human FER Protein
Catalog #: PEH7484



Summary

Name	FER_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P16591
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NNFIDVHRTT AAKEQEIEFD TSLLEENENL QANEIMWNNL TAESLQVMLK TLAEELMQTQ QMLLNKEEAV LELEKRIEES SETCEKKSDI VLLLSQKQAL EELKQSVQQL RCTEAKFSAQ KELLEQKVQE NDGKEPPPVV NYEEDARSVT SMERKERLSK

image



Background

Alternative Names	Tyrosine-protein kinase Fer (EC 2.7.10.2) (Feline encephalitis virus-related kinase FER) (Fujinami poultry sarcoma/Feline sarcoma-related protein Fer) (Proto-oncogene c-Fer) (Tyrosine kinase 3) (p94-Fer)
--------------------------	---

Note

Product Name: Recombinant Human FER Protein
Catalog #: PEH7484



For Research Use Only , Not for Diagnostic Use