


Product Name: Recombinant Human HSPA6 Protein
Catalog #: PEH7486



Summary

Name	HSP76_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P17066
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QDLLLDVAP LSLGLETAGG VMTTLIQRNA TIPTKQTQTF TTYS DNQPGV FIQVYEGERA MTKDNNLLGR FELSGIPPAP RGVPQIEVTF DIDANGILSV TATDRSTGKA NKITITNDKG RLSKEEVERM VHEAEQYKAE DEAQRDRVAA KNSLEAHVFH VKGSLQEESL RDKIPEEDRR KMQDKCREVL AWLEHNQLAE KEEYEHQKRE LEQICRPIFS RLYGGPGVPG GSSCGTQARQ GDPSTGPIIE

image



Background

Alternative Names Heat shock 70 kDa protein 6 (Heat shock 70 kDa protein B')

Note

Product Name: Recombinant Human HSPA6 Protein
Catalog #: PEH7486



For Research Use Only , Not for Diagnostic Use