


Product Name: Recombinant Human PFKL Protein
Catalog #: PEH7490



Summary

Name	PFKAL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P17858
Host	E.coli
Species	Human
Molecular Mass	104 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LTSGGDAQGM NAAVRVTRM GIYVGAKVFL IYEGY EGLVE GGENIKQANW LSVSNIIQLG GTIIGSARCK AFTTREGRR AAYNLVQHGI TNL CVIGGDG SLTGANIFRS EWGSLLEELV AEGKISSETA RTYSHLNIAG LVGSIDNDFC GTDMTIGTDS ALHRIMEVID AITTAQSHQ RTFVLEVMGR HCGYLALVSA LASGADWLFI PEAPPEDGWE NFM CERLGET RSRG SRLNII IIAEGAIDRN GKPISSSYVK DLVVQRLGFD TRVTVLGHVQ RGGTPSAFDR ILSSKMGMEA VMALLEATPD TPACVVTLSG NQSVRLPLME CVQMTKEVQK AMDDKRFDEA TQLRGGSFEN NWN IYKLLAH QKPPKEKSNF SLAILNVGAP AAGMNAAVRS AVRTGISHGH TVYVVDHGF E GLAKGQVQEV GWH DVAGWLG RGG SMLGTKR TLPKGQLESI VENIRIYGIH ALLVVG GFEA YEGVLQLVEA RGRYEELCIV MCVIPATISN NVP GTDFSLG SDTAVNAAME SCDRIKQSAS GTKRRVFIVE TMGGYCGYLA TVTGIAVGAD

image



Product Name: Recombinant Human PFKL Protein
Catalog #: PEH7490



Background

Alternative Names

ATP-dependent 6-phosphofructokinase, liver type (ATP-PFK) (PFK-L) (EC 2.7.1.11) (6-phosphofructokinase type B) (Phosphofructo-1-kinase isozyme B) (PFK-B) (Phosphohexokinase)

Note

For Research Use Only , Not for Diagnostic Use