


Product Name: Recombinant Human ALOX12 Protein
Catalog #: PEH7491



Summary

Name	LOX12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P18054
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GAWLFSGSYN RVQLWLVGTR GEAELELQLR PARGEEEFD HDVAEDLGLL QFVRLRKHHW LVDDAWFCDR ITVQGGPACA EVAFPCYRWV QGEDILSLPE

image



Background

Alternative Names	Polyunsaturated fatty acid lipoyxygenase ALOX12 (EC 1.13.11.-) (Arachidonate (12S)-lipoyxygenase) (12S-LOX) (12S-lipoyxygenase) (EC 1.13.11.31) (Arachidonate (15S)-lipoyxygenase) (EC 1.13.11.33) (Linoleate (13S)-lipoyxygenase) (Lipoxin synthase 12-LO) (EC 3.3.2.-) (Platelet-type lipoyxygenase 12)
--------------------------	---

Note

Product Name: Recombinant Human ALOX12 Protein
Catalog #: PEH7491



For Research Use Only , Not for Diagnostic Use