


Product Name: Recombinant Human TNFAIP3 Protein
Catalog #: PEH7502



Summary

Name	TNAP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P21580
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YLSNMRKAVK IRERTPEDIF KPTNGIIHHF KTMHRYTLEM FRTCQFCPQF REIIHKALID RNIQATLESQ KKLNWCREVR KLVALKTNGD GNCLMHATSQ YMWGVQDSDL VLRKALFSTL KETDTRNFKF RWQLESLSQ EFVETGLCYD TRNWNDEWDN LIKMASTDTP MARSGLQYNS LEEIHIFVLC NILRRPIIVI SDKMLRSLES GSNFAPLKV GGIYLPLHWP QECYRYPIVL GYDSHHFVPL VTLKDSGPEI RAVPLVNRDR GRFEDLKVHF LTDPENEMKE KLLKEYLMVI EIPVQGDHGT TTHLINAACL DEANLPKEIN LVDDYFELVQ HEYKKWQENS

image



Background

Alternative Names Tumor necrosis factor alpha-induced protein 3 (TNF alpha-induced protein 3) (EC 2.3.2.-) (EC 3.4.19.12) (OTU domain-containing protein 7C) (Putative)

Product Name: Recombinant Human TNFAIP3 Protein
Catalog #: PEH7502



DNA-binding protein A20) (Zinc finger protein A20) [Cleaved into: A20p50;
A20p37]

Note

For Research Use Only , Not for Diagnostic Use