


Summary

Name	UBA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P22314
Host	E.coli
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	<pre>RVSGPDPKPG SNCSPAQSVL SEVPSVPTNG MAKNGSEADI DEGLYSRQLY VLGHEAMKRL QTSSVLVSGL RGLGVEIAKN IILGGVKAVT LHDQGTQAWA DLSSQFYLRE EDIGKNRAEV SQPRLAELNS YVPVTAYTGP LVEDFLSGFQ VVVLTNTPLE DQLRVGEFCH NRGIKLVVAD TRGLFGQLFC DFG EEMILTD SNGEQPLSAM VSMVTKDNPG VVTCLDEARH GFESGDFVSF SEVQGMVELN GNQPMEIKVL GPYTFICDT SNFSDYIRGG IVSQVKVPPK ISFKSLVASL AEPDFVVTDF AKFSRPAQLH IGFQALHQFC AQHGRPPRPR NEEDAAELVA LAQAVNARAL PAVQQNNLDE DLIRKLAYVA AGDLAPINAF IGLLAAQEVN KACSGKFMPI MQWLYFDALE</pre>

image



Background

Product Name: Recombinant Human UBA1 Protein
Catalog #: PEH7507



Alternative Names

Ubiquitin-like modifier-activating enzyme 1 (EC 6.2.1.45) (Protein A1S9)
(Ubiquitin-activating enzyme E1)

Note

For Research Use Only , Not for Diagnostic Use