


Product Name: Recombinant Human SP100 Protein
Catalog #: PEH7512



Summary

Name	SP100_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P23497
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLEQGTGENS FRSLTWPPSG SPSHAGTTPP ENGLSEHPCE TEQINAKRKD TTSDKDDSLQ SQQTNEQCAQ KAEPTESCEQ IAVQVNNNGDA GREMPCPLPC DEESPEAELH NHGIQINSCS VRLVDIKKEK PFSNSKVECQ AQARTHNNQA SDIIVISSED SEGSTDVDEP LEVFISAPRS EPVINNDNPL ESNDEKEGQE ATCSRQIVP EPMDFRKLST FRESFKKRVV GQDHDFSESS EEEAPAEASS GALRSKHGEK APMTSRSTST WRIPSRKRRF SSSDFSDLN GEELQETCSS SLRRGSGSQP QEPENKKCSC VMCFPKGVPR SQEARTESSQ ASDMMDTMDV ENNSTLEKHS GKRRKKRRHR SKVNGLQRGR KKDRPRKHLT

image



Background

Alternative Names Nuclear autoantigen Sp-100 (Nuclear dot-associated Sp100 protein)

Product Name: Recombinant Human SP100 Protein
Catalog #: PEH7512



(Speckled 100 kDa)

Note

For Research Use Only , Not for Diagnostic Use