


Product Name: Recombinant Human NR2F2 Protein
Catalog #: PEH7515



Summary

Name	COT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P24468
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QTPGQGGPAS TPAQTAAGGQ GGPGGPGSDK QQQQHIIECV VCGDKSSGKH YGQFTCEGCK SFFKRSVRRN LSYTCRANRN CPIDQHHRNQ CQYCRLKKCL KVGMRREAVQ

image



Background

Alternative Names	COUP transcription factor 2 (COUP-TF2) (Apolipoprotein A-I regulatory protein 1) (ARP-1) (COUP transcription factor II) (COUP-TF II) (Nuclear receptor subfamily 2 group F member 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NR2F2 Protein
Catalog #: PEH7515

