



Product Name: Recombinant Human BRD2 Protein
Catalog #: PEH7516



Summary

Name	BRD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P25440
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLQDVFEFRY AKMPDEPLEP GPLPVSTAMP PGLAKSSSES SSEESSSESS SEEEEEDEEE DEEEEESESS DSEEEAHRL AELQEQLRAV HEQLAALSQG PISKPKRKRE KKEKKKKRKA EKHRGRAGAD EDDKGPRAPR PPQPKKSKKA SGSGGGSAAAL GPSGFGPSGG SGTKLPPKAT KTAPPALPTG YDSEEEEEER PMSYDEKRQL SLDINKLPGE KLGRVVHIIQ AREPSLRDSN PEEIEIDFET LKPSTLRELE RYVLSCLRKK PRKPYTIKKP VGKTKEELAL EKKRELEKRL QDVSGQLNST KKPPKKANEK TESSAQQVA VSRLSASSSS SDSSSSSSSS SSSDTSDSDS

image



Background

Alternative Names	Bromodomain-containing protein 2 (O27.1.1) (Really interesting new gene 3
--------------------------	---

Product Name: Recombinant Human BRD2 Protein
Catalog #: PEH7516



protein)

Note

For Research Use Only , Not for Diagnostic Use