


Product Name: Recombinant Human CD82 Protein
Catalog #: PEH7526



Summary

Name	CD82_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P27701
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VTELIRDYNS SREDSLQDAW DYVQAQVKCC GWVSFYNWTD NAELMNRPEV TYPSCSEVKG EEDNSLSVRK GFCEAPGNRT QSGNHPEDWP VYQEGCMEKV

image



Background

Alternative Names	CD82 antigen (C33 antigen) (IA4) (Inducible membrane protein R2) (Metastasis suppressor Kangai-1) (Suppressor of tumorigenicity 6 protein) (Tetraspanin-27) (Tspan-27) (CD antigen CD82)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use