


Product Name: Recombinant Human CAD Protein
Catalog #: PEH7528



Summary

Name	PYR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	P27708
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RPELGSRQDV EALWENMAVI DCFASDHAPH TLEEKCGSRP PPGFPGLETM LPLLLTAVSE GRLSLDDLLQ RLHHNPRRIF HLPPQEDTYV EVDLEHEWTI PSHMPFSKAH WTPFEGQKVK GTVRRVVLRG EVAYIDGQVL VPPGYGQDVR KWPQGAVPQL PPSAPATSEM TTTPERPRRG IPGLPDGRFH LPPRIHRASD PGLPAEEPKE KSSRKVAEPE LMGTPDGTCY PPPPVPRQAS

image



Background

Alternative Names	CAD protein [Includes: Glutamine-dependent carbamoyl-phosphate synthase (EC 6.3.5.5); Aspartate carbamoyltransferase (EC 2.1.3.2); Dihydroorotase (EC 3.5.2.3)]
--------------------------	---

Product Name: Recombinant Human CAD Protein
Catalog #: PEH7528



Note

For Research Use Only , Not for Diagnostic Use