


Summary

Name	PSB9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P28065
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TTIMAVEFDG GVVMGSDSRV SAGEAVVNRV FDKLSPLHER IYCALSGSAA DAQAVADMAA YQLELHGIEL EEPPLVLAAA NVVRNISYKY REDLSAHLMV AGWDQREGGQ VYGTLLGMLT RQPFAIGGSG STFIYGYVDA AYKPGMSPEE CRRFTTDAIA LAMSRDGSSG GVIYLVTTA AGVDHRVILG NELPKFYDE

image



Background

Alternative Names	Proteasome subunit beta type-9 (EC 3.4.25.1) (Low molecular mass protein 2) (Macropain chain 7) (Multicatalytic endopeptidase complex chain 7) (Proteasome chain 7) (Proteasome subunit beta-1i) (Really interesting new gene 12 protein)
--------------------------	---

Product Name: Recombinant Human PSMB9 Protein
Catalog #: PEH7529



Note

For Research Use Only , Not for Diagnostic Use