


Product Name: Recombinant Human TPP2 Protein
Catalog #: PEH7531



Summary

Name	TPP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P29144
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SEKEKDLKEE FTEALRDLKI QWMTKLDSSD IYNELKETYP NYLPLYVARL HQLDAEKERM KRLNEIVDAA NAVISHIDQT ALAVYIAMKT DPRPDAATIK NDMDKQKSTL VDALCRKGCA LADHLLHTQA QDGAISTDAE GKEEEGESPL DSLAETFWET TKWTDLFDNK VLTfAYKHAL

image



Background

Alternative Names	Tripeptidyl-peptidase 2 (TPP-2) (EC 3.4.14.10) (Tripeptidyl aminopeptidase) (Tripeptidyl-peptidase II) (TPP-II)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TPP2 Protein
Catalog #: PEH7531

