



Product Name: Recombinant Human CDC25A Protein
Catalog #: PEH7536



Summary

Name	MPIP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P30304
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PLDSKENLEN PMRRIHSLPQ KLLGCSPALK RSHSDSLDHD IFQLIDPDEN KENEAFFKK PVRPVSRLGCL HSHGLQEGKD LFTQRQNSAP ARMLSSNERD SSEPGNFIPL FTPQSPVTAT LSDEDDGFVD LLDGENLKNE EETPSCMASL WTAPLVMRTT NLDNRCKLFD SPSLCSSSTR SVLKRPERSQ EESPPGSTKR RKSMMSGASPK ESTNPEKAHE TLHQSLSLAS

image



Background

Alternative Names	M-phase inducer phosphatase 1 (EC 3.1.3.48) (Dual specificity phosphatase Cdc25A)
--------------------------	---

Note

Product Name: Recombinant Human CDC25A Protein
Catalog #: PEH7536



For Research Use Only , Not for Diagnostic Use